



Better access^{†,‡}

than the Harmonic Focus^{TM*} + because of its optimized form factor and fine, curved jaw profile

More effective blunt dissection^{†,‡}

than the Harmonic Focus^{TM*} + due to its symmetric jaw, which applies equal force on opposing tissue planes

THE COOLER, MORE PRECISE DISSECTOR.

The LigaSureTM exact dissector outperforms the Harmonic Focus^{TM*} + in several key areas

It's All About the Jaw

†11 out of 12 surgeons surveyed agreed.
‡11 out of 13 surgeons surveyed agreed.


Medtronic
Further, Together

Cooler Jaw Temperature²

The nonstick nano-coated jaws on the LigaSure™ exact dissector are:

177°C

cooler than the Harmonic Focus™+ after 5 activations



Jaw Temperature After One Activation²


<60°C

LigaSure™ exact dissector

221°C

Harmonic Focus™+

Faster Cool Down Times²



LigaSure™ Exact Dissector

0 seconds

to cool down to 60°C after 1 activation

10 seconds

to cool down to 60°C after 5 activations



Harmonic Focus™+

32 seconds

to cool down to 60°C after 1 activation

36 seconds

to cool down to 60°C after 5 activations

BETTER ACCESS.¹ MORE PRECISE DISSECTION.¹ FASTER COOLING TIMES.²

The LigaSure™ exact dissector provides all three. When compared to the Harmonic Focus™+, it delivers:

- Fewer instrument exchanges[‡]
- A more efficient procedure[†]
- The reliable performance of the Valleylab™ FT10 and LS10 (where available) energy platforms

[†] 11 out of 12 surgeons surveyed agreed.

[‡] 11 out of 13 surgeons surveyed agreed.

[‡] 13 out of 14 surgeons surveyed agreed.

1. Based on internal report #RE00114823. Validation labs: surgeon evaluation of LigaSure™ exact dissector, nano-coated. July 18 to 26, 2017.

2. Based on internal report #RE00107711 Rev A. Market research: Thermal profile comparison of LigaSure™ exact dissector and the Ethicon Harmonic Focus™+ conducted on porcine tissue. Aug. 29, 2017.

© 2020 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further. Together are trademarks of Medtronic. ™* Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company. 04/2020-US171002a(1)–[WF# 2322960]

Let's bring a cooler,² more precise dissection^{1,2} device to your OR.

For more information, call your local sales representative or visit [medtronic.com/covidien](https://www.medtronic.com/covidien)

Medtronic